

# Give Patients Every Opportunity For Success

## Luna G3, React Health Masks and V-Com™

The Luna G3 PAP and Bilevel Devices are designed for patient compliance with patient forward standard features like Full Color LED Screen, Standard Heated Humidifier, Auto Ramp Feature, Supply Reminders and Bedside clock.



The Siesta Full Face, Siesta Nasal and Rio II Nasal Pillows Mask have a less than 1% refit rate\*

Comfort-related complaints are among the top reasons for noncompliance to PAP therapy.\* The V-Com is a flow attenuator placed in between the mask and CPAP device, designed to provide comfort for patients acclimating to PAP therapy.

\* Aggregate independent data from SleepGlad Mask Filter May 202-April 2021

\* Bollig SM. Encouraging CPAP adherence: it is everyone's job. Respir Care. 2010 Sep;55(9):1230-9. PMID: 20800003.

## Luna G3 PAP Devices

The Luna G3 CPAP (LG3500) APAP (LG3600) and the Luna G3 Bilevel 25A have all the standard features you have come to expect in a device in a smaller, more modern footprint.



- » Luna G3 Devices include:
- » Full Color LED Screen
- » Standard Heated Tubing
- » Ramp time (0-60 mins)
- » Cellular Modems Available on all devices
- » Bedside Clock
- » Supply Reminders
- » And more!



## React Health Masks

### Siesta Full Face Mask

- » No forehead brace allows for increased line of sight - even with your glasses on
- » Self-adjusting ergonomic cushion design allows for superior seal
- » Padded support arms minimize irritation and marks around cheekbones and eyes
- » Comes in a fit pack and individual size configuration
- » Component Parts also available separately



### Siesta Nasal Mask

- » Padded support arms minimize irritation and marks around cheekbones and eyes Quick release clips allow for simple attachment and removal of mask
- » Easy change cushion cup snaps into place allows for ease of removal and cleaning
- » Comes in a fit pack and individual size configurations
- » Component Parts also available separately



### Rio II Nasal Pillows Mask

- » A lightweight experience- Just 2.3 ounces: One of the lightest masks on the market today
- » Dual wall pillows with trampoline bounce for secure fit
- » Split back strap allows for a better fit, while
- » Velcro fasteners create added layer of adjustment
- » Comes in a convenient fit pack including S,M, L Cushions
- » Component Parts also available separately



## V-Com™ Comfort Accessory

Compliance to therapy is one of the most important factors when starting a patient on PAP therapy and can also be a clinician's greatest challenge. Patients often complain that PAP therapy is not comfortable. Patients who are less comfortable may not gain early acceptance to therapy, which has been shown to increase long term adherence\* Adding the V-Com to your Luna G3 PAP and React Health mask can be just what the patient needs to get acclimated to their PAP therapy!\*\*

The V-Com is a flow attenuator placed in between the mask and PAP device, designed to provide comfort for patients acclimating to CPAP therapy.

V-com™ softens peak inspiratory flow which often overwhelms new PAP users and provides inspiratory comfort with minimal effect on expiratory positive airway pressure (EPAP) and therapy.

*\* Budhiraja R et al. Early CPAP use identifies subsequent adherence to CPAP therapy. SLEEP 2007;30(3):320-324*

*\*\* Consult with the package insert to determine patient eligibility for V-Com*

